



1.5 ml, Screw Neck Vial N8, 32 x 11.6 mm,
clear glass, 1st hydrol. class,
small opening, label + filling lines



1.5 ml, Screw Neck Vial N8, 32 x 11.6 mm,
amber glass, 1st hydrol. class,
small opening, label + filling lines



8 mm Black Screw Cap with Septa
White PTFE Red Silicone



1.5 ml, PP Short Thread Vial, transparent,
with filling lines, 32 x 11.6 mm,
slightly concave shaped bottom



1.5 ml, Short Thread Vial ND9, 32 x 11.6 mm,
clear glass, 1st hydrol. class,
wide opening, label + filling lines



1.5 ml, Short Thread Vial ND9, 32 x 11.6 mm,
amber glass, 1st hydrol. class,
wide opening, label + filling lines



9 mm Blue Screw Cap with Septa
White PTFE Red Silicone



9 mm PP Short Thread Cap Blue,
6 mm centre hole (Silicone
beige / PTFE with slit bonding)



0.3 ml PP Short Thread Micro Vial,
transparent, 32 x 11.6 mm



1.5 ml, Crimp Neck Vial, 32 x 11.6 mm,
clear glass, 1st hydrol. class,
wide opening, label + filling lines



1.5 ml, Crimp Neck Vial, 32 x 11.6 mm,
amber glass, 1st hydrol. class,
wide opening, label + filling lines



11 mm Aluminium Crimp Cap
Silicone white / PTFE red Septa



20ml Headspace-Vial,
75.5 x 22.5mm, clear glass,
1st hydrol. class, DIN-crimp neck,
long neck, flat bottom



20ml Headspace-Vial,
75.5 x 22.5mm, amber glass,
1st hydrol. class, DIN-crimp neck,
long neck, rounded bottom



10 ml Headspace Vial,
46 x 22.5 mm clear glass,
1st hydrol. class, DIN crimp neck,
rounded bottom



10 ml Headspace Vial,
46 x 22.5 mm amber glass,
1st hydrol. class, DIN crimp neck,
long neck, flat bottom



Conical Insert with
Plastic Spring



Flat Bottom
Insert



20 mm Aluminium
Crimp Cap with PTFE
Silicon Septa



20 mm Magnetic
Crimp Cap with
PTFE Silicon Septa



20 mm Bimetallic
Crimp Cap with
PTFE Silicon Septa



20 ml Headspace Vial,
Clear Glass, Screw Neck



18 mm Magnetic Screw with Cap
(8 mm centre hole) with septa



40 ml Clear Vial with
Screw Cap and Septa



1 ml Shell Vial, 40 x 8.2 mm, clear glass, 1st hydro. class,
8 mm PE-Plug, soft, without insertion barrier for
Micro-insert, transparent for Waters
Wisp 96 Pos. Carousel, Shimadzu



20 ml Vial Rack



2 ml Vial Rack - 50 slots



2 ml Vial Rack
with Lid - 100 slots

Vial Rack



10 ml Glass Vials, Amber Colored, Crimp Top,
Rubber Stopper with
20 mm Aluminium Seal Cap



10 ml Amber Colour Storage Vial with
Black Screw Cap Closed with
Inner Plug



20 ml Glass Vials, Amber Colored, Crimp Top,
Rubber Stopper with
20 mm Aluminium Seal Cap