

## GC& HPLC SYRINGES

- Fixes Needle Syringes for GC
- Removal Needle Syringes for GC
- Fixed Needle Syringes Flexible Plunger
- Removable Syringes Flexible Plunger
- Guided Plunger Fixed Needle Syringes
- Guided Plunger Removable Needle Syringes
- Gas Tight Syringes Fixed Needle for GC
- Gas Tight Syringes Removable Needle for GC
- SGE Repeating Adaptor Syringes
- Luer-Lock Gas-Tight Syringes
- For Rheodyne and Valco Valve - HPLC
- Manual Gas Syringes - Plunger in Needle Syringes
- SGE Removable Needle



## GC Capillary Columns

- BP-1 100% Dimethyl Polysiloxane
- BP-5 5% Phenyl Polysiloxane
- BP-10 (1701) - 14% Cyanopropylphenyl / Polysiloxane
- BP-20 (WAX) Polyethylene Glycol
- BP-21 (FFAP) Polyethylene Glycol (TPA Treated)
- BPX-5 5% Phenyl (Equiv.) Polysilphenylene Siloxane
- BP-225 - 50% Cyanopropylphenyl Polysiloxane
- BP-1 Pona 100% Dimethyl Polysiloxane
- Cydex-B Permethyl Ateed Beta Cyclodextrin
- BP-624 Cyanopropylphenyl Polysiloxane
- BPX-35 - 35% Phenyl (Equiv.) Polysilphenylene - Siloxane
- BPX-50 - 50% Phenyl (Equiv.) Polysilphenylene - Siloxane
- BPX-70 - 70% Cyanopropyl (Equiv.) Polysilphenylene - Siloxane



## GC Graphite Ferrule

- 1/16" Fittings - All Ferrules supplied in 100% Graphite
- 1/16" Fittings - All Ferrules supplied in 15% Graphite / 85% VESPEL



## Inlet Liners for Agilent, Thermo Instruments



Split / Splitless FocusLiner



Split/ Splitless with single taper (quartz wool)



Split/ Splitless quartz, straight-through liner

## Inlet Liners for PerkinElmer Instruments



Split / Splitless FocusLiner



Split / Splitless tapered FocusLiner



Split, straight-through liner



Split / Splitless for PSS injector



Split / Splitless FocusLiner for PSS injector

## Inlet Liners for Shimadzu Instruments



Split / Splitless FocusLiner



Split/ Splitless with middle gooseneck



Split/ Splitless with middle gooseneck



Connectite (0.53 mm ID columns)

## GC Septa

GP Grade	EC Grade	MN Grade	HT Grade	Enduro Blue