

## GC COLUMN SAFETY CABINETS

# PcolSafe®-GC

(Column Safety Cabinet)

- Safe column organizer cabinet for capillary columns
- 5 Drawers for each capillary column
- Each drawer provided with soft cushioning for safety of capillary columns from jerks & shocks
- MOC : MS Powder coated
- Central Lock

### Specification

#### ■ CSC 5

- Capacity : 5 GC column
- Dimensions (WxHxD) : 295 x 590 x 410 mm
- No. of Drawers : 5 drawers

#### ■ CSC 10

- Capacity : 10 GC column
- Dimensions (WxHxD) : 850 x 600 x 420 mm
- No. of Drawers : 5 drawers (2 Capillary columns storage per drawer)



MODEL  
CSC 5



MODEL  
CSC 10