

NITROGEN SAMPLE EVAPORATOR / CONCENTRATOR WELL PLATE



Model : EV-PLUS 96 X 2 WP

- * **Capacity:** Two Standard 96 x 2 well plates
- * **Max. sample volume:** 1 mL (as per WP specs)
- * **Well plate Bay:** 2 included as per requirement
- * **Number of Samples:** Single Well plate of 96 samples X 2
- * **Heating:** Block Heating for each plate
- * **Timer Range:** 1 to 99 minutes/indefinite Gas Regulator and Gas Gauge Range: 0 to 20 psi
- * **Gas Supply Requirements:** Min. inlet pressure 60 psi, Max.inlet pressure 100 psi
- * **Exhaust Port:** In-built , Duct hose provided
- * **Power Requirements:** 230 VAC at 50 Hz
- * **Technology:** Vortex technology
- * **Gas control:** Single input
- * **Net Weight:** 15 kg Approx